

## Antifade

Section of Cancer Genomics, Genetics Branch, NCI  
National Institutes of Health

### Antifade 1: 1,4-phenylene-diamine (This is the preferred antifade reagent)

#### Reagents

**1,4-phenylene-diamine**  
Sigma, Cat. P1519, 100 g  
**Glycerol, 86%**  
**Phosphate Buffered Saline (PBS), 1X**  
**Sodium bicarbonate**  
**Sodium carbonate**

#### Preparation

**Carbonate-Bicarbonate Buffer (pH 9.0)**  
Sodium bicarbonate, 0.5 M (pH 8.13)      4 ml  
Sodium carbonate, 0.5 M (pH 11.32)      1 ml  
Filter sterilize

#### Procedure

1. Prepare carbonate-bicarbonate buffer (pH 9.0).
2. Dissolve 50 mg 1,4-phenylenediamine in 2 ml 1X PBS.
3. Adjust pH with carbonate-bicarbonate buffer to 8.0.
4. Add 1X PBS to 5 ml. **Important:** Add 2 ml buffer and check pH, then add dropwise until pH is 7.99-8.00. If pH exceeds 8.0, the procedure must be started over.
5. Mix with 45 ml 86% glycerol. Leave on inverter for at least 1 hr.
6. Aliquot and store at -20°C.

**Antifade 2: DABCO (1,4 diazabicyclo[2.2.2]octane)****Reagents****DABCO (1,4-diazabicyclo[2.2.2]octane)**

Sigma, Cat. D2522

**Glycerol, 86%****Tris-HCl, 1 M pH 8.0****Water, sterile****Preparation**

<b><u>Components</u></b>	<b><u>Amount</u></b>
DABCO	0.233 g
Tris-HCl, 1 M pH 8.0	200 $\mu$ l
Sterile water	800 $\mu$ l
Glycerin, 86%	9 ml

**Procedure**

1. Combine components.
2. Dissolve by warming to 70°C.
3. Vortex.
4. Aliquot and store at -20°C.