

cDNA Synthesis for RT-PCR Protocol

Section of Cancer Genomics, Genetics Branch, NCI
National Institutes of Health

Reagents

5X First Strand Buffer

Invitrogen, Cat. Y00146

10mM dNTP Set

Amersham Biosciences, Cat. 27-2035-01

0.1M DTT

Invitrogen, Cat. Y00147

Random Primers

Invitrogen, Cat. 48190-011

RNaseOUT Ribonuclease Inhibitor

Invitrogen, Cat. 10777-019

SuperScript II RNase H- Reverse Transcriptase

Invitrogen, Cat. 18064-014

Procedure

1. Combine 3-5 μg total RNA and molecular grade water to 8 μl final volume.
2. Add 3 μl Random Primers.
3. Add 1 μl dNTP mix.
4. Vortex and then spin down tube.
5. Incubate at 65°C for 5 min.
6. Place tube on ice.
7. Add 4 μl 5X Buffer, 2 μl DTT and 1 μl RNaseOut.
8. Vortex and then spin down tube.
9. Incubate at 42°C for 1 min.
10. Add 1 μl SuperScript II RNase H- Reverse Transcriptase.
11. Incubate at 42°C for 60 min.
12. Incubate at 70°C for 15 min.
13. Add 180 μl molecular grade water.

14. Use Nanodrop 1000 to measure concentration. Set sample type setting to Other Sample and the constant to 33.
15. Store at -80°C.

NOTES

1. 10 mM dNTP Mix:
 - 10 μ l dATP [100mM]
 - 10 μ l dCTP [100mM]
 - 10 μ l dGTP [100mM]
 - 10 μ l dTTP [100mM]
 - 60 μ l molecular grade water