

DNA Precipitation and Hybridization (CGH)

Section of Cancer Genomics, Genetics Branch, NCI
National Institutes of Health

Reagents

Dextran sulfate (50%)

Intergen, Cat. S4030

Ethanol, absolute

Formamide, deionized

Mouse Cot-1TM DNA, Cat. 18440-016, 500 mg

Human Cot-1TM DNA, 1 mg/ml

Invitrogen Corp., Cat. 15279-011, 500 µg

Salmon testes DNA, 9.7 mg/ml

SIGMA Molecular Biology, Cat. D-7657, 1 ml

SSC, 20X

Sodium acetate (Na-Acetate), 3M

Water, sterile

Preparation

Master Mix

Dextran sulfate, 50%	40 ml
----------------------	-------

20X SSC	20 ml
---------	-------

Sterile H ₂ O	40 ml
--------------------------	-------

Total	100 ml
--------------	---------------

Vortex solution and place tube on a shaking platform overnight to insure proper mixing.

Aliquot, and store at -20°C

70% Formamide/2X SSC

Deionized formamide	70 µl
---------------------	-------

20X SSC	10 µl
---------	-------

Sterile water	20 µl
---------------	-------

Total	100 µl
--------------	---------------

Adjust to pH 7.0

Procedure

1. Add to an eppendorf tube:
 - 10-25 μ l nick-translated test or tumor probe DNA (500-1000 ng DNA)
 - Equal amount of nick-translated control whole genomic DNA as probe DNA (Note: can use 1-2 μ g DNA if tumor DNA is isolated from paraffin material)
 - 30-60 μ l human Cot-1 DNA (1 mg/ml)
 - 1 μ l salmon sperm DNA (10 mg/ml)
- Note:** Usually the test DNA is nick translated with Biotin-16-dUTP and the control DNA is nick translated with Digoxigenin-11-dUTP.
2. Add 1/10 volume Na-Acetate (3M).
3. Add 2.5-3.0 x total volume of absolute ethanol.
4. Vortex, store at -20°C overnight, or at -80°C for at least 15-30 min.
5. Centrifuge (14000 rpm) precipitated DNA at 4°C for 30 min.
6. Pour off supernatant and speed vac for 5-10 min to dry pellet.
7. Add 5 μ l pre-warmed deionized formamide (pH 7.5), incubate at 37°C for 30 min, shaking; vortex a few times during the 30 min incubation.
8. Add 5 μ l pre-warmed Master Mix, vortex, and centrifuge briefly. Incubate at 37°C 15-30 min.
9. Denature probe DNA at 80°C for 5 min and centrifuge briefly. Probe can be kept at 37°C until ready to denature.
10. Preanneal at 37°C for 1-2 hr.
11. For slide denaturation apply 120 μ l 70% formamide/2X SSC to a 24 x 60 mm coverslip and place inverted slide onto coverslip.
12. Incubate slides at 75°C for 1.5 min on a slide warmer.
13. Remove coverslip and immediately place slide in ice cold 70% ethanol for 3 min, followed by 90% ethanol and 100% ethanol for 3 min each; air dry.
14. Add pre-annealed probe DNA to denaturated slides and cover with 18 mm² or 22 mm² coverslips; seal coverslips with rubber cement.
15. Hybridize at 37°C in a humidified chamber for 72 hr.