

High Resolution Human Lymphocyte Chromosomes

Section of Cancer Genomics, Genetics Branch, NCI
National Institutes of Health

Reagents

Acetic acid, glacial

Mallinckrodt, Cat. V193

Bromodeoxyuridine (BrdU)

Sigma, Cat. B9285

Colcemid, KaryoMAX Colcemid Solution (10 µg/ml)

Invitrogen Corp., Cat. 15210-016

Fetal Bovine Serum (FBS) Qualified, heat inactivated

Invitrogen Corp., Cat. 16140-022, 500 ml

L-Glutamine-200 mM, 100x

Invitrogen Corp., Cat. 25030-016

Penicillin/Streptomycin 5,000 U/ml/5,000 µg/mL

Invitrogen Corp., Cat. 15070-014

Phytohaemagglutinin (PHA), HA 15

Murex Diagnostics Ltd., Dartford, England DA1 5LR

Methotrexate (MTX)

Sigma, Cat. M8407

Methyl alcohol, anhydrous

Mallinckrodt, Cat. 3016

Potassium chloride (KCl)

Mallinckrodt, Cat. 6858

RPMI Medium 1640

Invitrogen Corp., Cat. 21870-050

Preparation

RPMI 1640 Complete Medium

Components	Amount
RPMI Medium 1640	385 ml
L-Glutamine-200 mM, 100x	5 ml
Penicillin/Streptomycin 5,000 U/ml/5,000 µg/ml	10 ml
Fetal Bovine Serum Qualified, heat inactivated	100 ml

Hypotonic solution: 0.075M KCl

KCl	5.6 g
Distilled water	1000 ml

Fixative

Methanol/glacial acetic acid, 3:1 (volume:volume)

MTX stock (prepare fresh)

10^{-5} M in H₂O

BrdU stock (light sensitive)

1 mg/ml in distilled water

Procedure

1. Use T25 (with 5 ml media) or T75 flasks (with 20 ml media).
2. Initiate PHA-stimulated lymphocyte cultures. Incubate in upright position at 37°C.
3. At 72 h add from the MTX stock (10^{-5} M) to a final concentration of 10^{-7} M; mix well and incubate an additional 17 hr.
4. After 17 h centrifuge the contents of the flasks, remove the supernatant, and wash the pellet twice with unsupplemented media.
5. After the second wash resuspend the pellet in RPMI 1640 20% FBS and transfer to a fresh flask.
6. Add from the BrdU stock (1 mg/ml) to a final concentration of 25 µg/ml (minimize light exposure).
7. Incubate for 5 h 30 min at 37°C.
8. During the last 10 min. of incubation add Colcemid stock (10 µg/ml) to a final concentration of 0.06 µg/ml.
9. Centrifuge cultures for 10 min.
10. Remove supernatant and add hypotonic solution, 0.075M KCl; incubate for 10 min at 37°C.
11. Add a few drops of fresh fixative and spin for 10 min at 1,000 rpm.
12. Aspirate supernatant and add fixative.
13. Resuspend pellet gently in fixative and centrifuge for 5 min at 1,200 rpm. Repeat this step two more times.

14. Store pellet under fixative at -20°C until ready to prepare slides.