

Nick Translation (CGH)

Section of Cancer Genomics, Genetics Branch, NCI
National Institutes of Health

Reagents

Biotin-16-dUTP

Boehringer Mannheim, Cat. 1093070

Bovine serum albumin (BSA)

Dig-11-dUTP

Boehringer Mannheim, Cat. 1093088

dATP, dTTP, dGTP, dCTP

Boehringer Mannheim, Cat. 105 1440, 105 1458, 105 1466, 105 1482

DNase I from bovine pancreas

Boehringer Mannheim, Cat.104 159, 100 mg

EDTA, 0.5 M

Glycerol

Lambda HindIII DNA marker

Magnesium chloride (MgCl₂), 2 M

β-Mercaptoethanol, 99%

Polymerase (Kornberg)

Boehringer Mannheim, Cat. 104 485

NaCl, 1 M

Tris-HCl, 1 M, pH 8.0

Water, sterile

Preparation

		<u>final conc.</u>
dNTP		
100 mM dATP, dCTP, and dGTP	5 µl of each	0.5 mM
100 mM dTTP	1 µl	0.1 mM
Sterile water	984 µl	
Total	1000 µl	
(equals 0.5 mM each of dATP, dCTP, and dGTP, and 0.05 mM dTTP)		
*Aliquot and store at -20°C		

		<u>final conc.</u>
DNase I stock solution, 1mg/ml		
DNase I	10 mg	1 mg/ml
NaCl, 1M	1.5 ml	0.15 M
Glycerol	5 ml	50%
Sterile water	bring up to 10 ml	
Total	10 ml	
*Aliquot and store at -20°C		

		<u>final conc.</u>
10X NT-Buffer		
Tris-HCL, 1 M, pH 8.0	500 µl	0.5 M

MgCl ₂ , 2 M	25 μ l	50 mM
BSA, 10 mg/ml	50 μ l	0.5 mg/ml
Sterile water	425 μ l	
Total	1 ml	

*Aliquot and store at -20°C

0.1M β-Mercaptoethanol		<u>final conc.</u>
99% solution (14.4 M)	34.7 μ l	0.1 M
Sterile water	bring up to 5 ml	
Total	5 ml	

*Aliquot and store at -20°C

Procedure

- For each DNA sample, add to an eppendorf tube:
 - 2 μ g DNA
 - 10 μ l 10X NT-Buffer
 - 10 μ l dNTP mix
 - 10 μ l 0.1 M β -Mercaptoethanol
 - 4 μ l BIO-16-dUTP or 4 μ l DIG-11-dUTP (1 mM)
 - X μ l sterile water
 (The total volume including reagents added in step 3 should be 100 μ l)
- Vortex, centrifuge, and place tubes on ice.
- Add 2 μ l Polymerase (Kornberg) first, and then 3-8 μ l Dnase (1 mg/ml) 1:1000 (Dnase amount will vary from different batches).
- Vortex and centrifuge and centrifuge.
- Incubate at 15°C for 2 hr (1.5-2 hr).
- Prepare gel electrophoresis.
- Run about 5 μ l of each sample and the Lambda HindIII DNA marker; ideally the length of the DNA should be 500-900 bp after nick translation.
- If DNA is too large, add more DNase and incubate at 15°C for 10-30 min.
- Stop the nick translation with 1 μ l of 0.5 M EDTA and incubate at 65°C for 10 min.
- Store DNA at -20°C or precipitate the same day.