

Nick Translation (FISH)

Section of Cancer Genomics, Genetics Branch, NCI
National Institutes of Health

Reagents

Bovine serum albumin (BSA)

dATP, dTTP, dGTP, dCTP

Boehringer Mannheim, Cat. 105 1440, 105 1458, 105 1466, 105 1482

DNase I from bovine pancreas

Boehringer Mannheim, Cat. 104 159, 100 mg

dUTP (conjugated to hapten or fluorochrome of choice)

EDTA, 0.5 M

Glycerol

Lambda HindIII DNA marker

Magnesium chloride (MgCl₂), 0.5 M

β-Mercaptoethanol, 99%

Polymerase (Kornberg)

Boehringer Mannheim, Cat. 104 485

NaCl, 1 M

Tris-HCl, 1 M, pH 8.0

Water, sterile

Preparation

dNTP

100 mM dATP, dCTP, and dGTP	5 μl of each	f.c. [0.5 mM]
-----------------------------	--------------	---------------

100 mM dTTP	1 μl	f.c. [0.1 mM]
-------------	------	---------------

Sterile water	984 μl	
---------------	--------	--

Total	1000 μl	
-------	---------	--

(equals 0.5 mM each of dATP, dCTP, and dGTP, and 0.05 mM dTTP)

*Aliquot and store at -20°C

DNase I stock solution, 1mg/ml

DNase I	10 mg	
---------	-------	--

NaCl, 1M	1.5 ml	f.c. [0.15 M]
----------	--------	---------------

Glycerol	5 ml	f.c. [50%]
----------	------	------------

Sterile water	bring up to 10 ml	
---------------	-------------------	--

*Aliquot and store at -20°C

10X NT-Buffer

Tris-HCL, 1 M, pH 8.0	500 μ l f.c. [0.5 M]
MgCl ₂ , 0.5 M	100 μ l f.c. [50 mM]
BSA, 10 mg/ml	50 μ l f.c. [0.5 mg.ml]
Sterile water	350 μ l
<hr/> Total	<hr/> 1000 μ l

*Aliquot and store at -20°C

0.1M β -Mercaptoethanol

99% solution (14.4 M)	34.7 μ l
Sterile water	bring up to 5 ml

*Aliquot and store at -20°C

Procedure

- For each DNA sample, add to an eppendorph tube:
 - 2 μ g DNA
 - 10 μ l 10X NT-Buffer
 - 10 μ l dNTP
 - 10 μ l 0.1 M β -Mercaptoethanol
 - 4 μ l BIO-16-dUTP or 4 μ l DIG-11-dUTP (1 mM)
 - X μ l sterile water
 (The total volume including reagents added in step 3 should be 100 μ l)
- Vortex, centrifuge, and place tubes on ice.
- Add 2 μ l Polymerase (Kornberg) first, and then 3-8 μ l Dnase (1 mg/ml) diluted 1:1000
- Flick tube to mix.
- Incubate at 15°C for 1.5-2 hr.
- Prepare gel electrophoresis.
- Run 5 μ l of each sample with loading buffer and the Lambda HindIII DNA marker; ideally the length of the DNA should be 500-900 bp for chromosome paint probes or 300-600 bp for gene specific probes after nick translation.
- If DNA is too large, add more DNase and incubate at 15°C for 10-30 min.
- Stop the nick translation with 1 μ l of 0.5 M EDTA and incubate at 65°C for 10 min.
- Store DNA at -20°C or precipitate the same day for hybridization.